

Utilities for Mass Spectrometry Analysis of Proteins

User's Manual

Version 2.2.0

May 2022

To download Utilities for Mass Spectrometry Analysis of Proteins visit:

www.umsap.nl

Utilities for Mass Spectrometry Analysis of Proteins
Copyright © 2017 Kenny Bravo Rodriguez.
All Rights Reserved.

Contents

List of Figures	III
List of Tables	IV
1 Introduction	1
1.1 Citing Utilities for Mass Spectrometry Analysis of Proteins	1
1.2 Acknowledgments	2
2 Obtaining and Installing Utilities for Mass Spectrometry Analysis of Proteins	3
2.1 Obtaining Utilities for Mass Spectrometry Analysis of Proteins	3
2.2 Installing Utilities for Mass Spectrometry Analysis of Proteins	3
2.3 Uninstalling Utilities for Mass Spectrometry Analysis of Proteins	4

List of Figures

List of Tables

Chapter 1

Introduction

Utilities for Mass Spectrometry Analysis of Proteins (UMSAP) is a graphical user interface (GUI) designed to speed up the post-processing of data obtained during mass spectrometry studies involving proteins. The program is not intended to analyze a mass spectrum or a mass chromatogram, neither to identify the peaks in a mass spectrum. The main objective is the fast post-processing of the vast amount of data generated in mass spectrometry experiments involving proteins after peak identification have been performed.

The program is organized in modules with each module performing a single type of data post-processing. The reason for this clear separation is the high dependency between the type of mass spectrometry experiment performed and the way in which the resulting data must be post-processed. The modules are designed in such a way that the required user input is minimized but still users can control every aspect of the analysis. Currently, the software contains three modules, but several others are already planned.

1.1 Citing Utilities for Mass Spectrometry Analysis of Proteins

If results obtained with UMSAP are published in any way, please acknowledge the use of UMSAP by including the following sentence:

”Utilities for Mass Spectrometry Analysis of Proteins was created by Kenny Bravo Rodriguez at the University of Duisburg-Essen and is currently developed at the Max Planck Institute of Molecular Physiology.”

Any published work, which uses UMSAP, should include the following reference:

Kenny Bravo-Rodriguez, Birte Hagemeier, Lea Drescher, Marian Lorenz, Michael Meltzer, Farnusch Kaschani, Markus Kaiser and Michael Ehrmann. (2018). Utilities for Mass Spectrometry Analysis of Proteins (UMSAP): Fast post-processing of mass spectrometry data. [Rapid Communications in Mass Spectrometry](#), 32(19), 1659–1667.

Electronic documents should include a direct link to the official web page of UMSAP at: www.umsap.nl

1.2 Acknowledgments

I would like to thank all the persons that have contributed to the development of UMSAP, either by contributing ideas and suggestions or by testing the code. Special thanks go to: Dr. Farnusch Kaschani, Dr. Juliana Rey, Dr. Petra Janning and Prof. Dr. Daniel Hoffmann.

In particular, I would like to thank Prof. Dr. Michael Ehrmann.

Chapter 2

Obtaining and Installing Utilities for Mass Spectrometry Analysis of Proteins

2.1 Obtaining Utilities for Mass Spectrometry Analysis of Proteins

UMSAP is distributed free of charge for anyone interested in using it. To obtain a copy of the software just register at www.umsap.nl and go to the Download page.

No extra software or packages are needed for UMSAP to properly work.

So far, UMSAP have been tested in macOS 10.12.6 and 10.14.6 and Windows 10. Linux users may download the source code of the software and adapt it to their specific distribution of Linux. Support for some Linux distributions will be available in the future.

2.2 Installing Utilities for Mass Spectrometry Analysis of Proteins

Windows

Unzip the file you just downloaded from www.umsap.nl. Then, copy the folder UMSAP to the location in your file system where you want to keep it. Finally, create a shortcut to the executable file UMSAP.exe found inside the main folder UMSAP. That is all. You are now ready to use UMSAP.

macOS

Unzip the file you just downloaded from www.umsap.nl. Then, just move the UMSAP.app folder to /Applications/. That is all. You are now ready to use UMSAP.

Depending on the security settings in macOS, it may be needed to explicitly allow UMSAP to be opened the first time the app is used.

2.3 Uninstalling Utilities for Mass Spectrometry Analysis of Proteins

UMSAP will not create any installation file in your computer. Therefore, the only thing you need to do, to completely uninstall UMSAP, is to delete the folder UMSAP.app in macOS or UMSAP in Windows. You should also delete any shortcut pointing to the executable file of UMSAP and the configuration file `.umsap_config.json` in your home folder. That is all.